Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

## 1 Identification of the substance/mixture and of the company/undertaking

- · Product identifier
- · Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)
- · Article number: 490-001.01R01, 4980, 13866
- · Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- · Application of the substance / the mixture Void
- Details of the supplier of the safety data sheet
- Manufacturer/Supplier:

Manufactured by NuSmile Ltd for:

Ultradent Products Inc.

505 W. Ultradent Drive (10200 S)

South Jordan, UT 84095-3942

USA

onlineordersupport@ultradent.com

EC Responsible Person

Ultradent Products GmbH

Am Westhover Berg 30

51149 Cologne Germany

 ${\it Email: infoDE@ultradent.com}$ 

Emergency Phone: +49(0)2203-35-92-0

- · Further information obtainable from: Customer Service
- · Emergency telephone number:

CHEMTREC (NORTH AMERICA) : (800) 424-9300 (INTERNATIONAL) : +(703) 527-3887

#### 2 Hazards identification

- · Classification of the substance or mixture
- · Classification according to Regulation (EC) No 1272/2008



Skin Corr. 1A H314 Causes severe skin burns and eye damage.

- · Additional information: Emergency Overview: Slightly acidic water-based gel solution
- · Label elements
- · Labelling according to Regulation (EC) No 1272/2008

The product is classified and labelled according to the GB CLP regulation.

- · Hazard pictograms GHS05
- · Signal word Danger
- · Hazard statements

H314 Causes severe skin burns and eye damage.

· Precautionary statements

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

*P103* Read carefully and follow all instructions.

P260 Do not breathe dusts or mists.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water

[or shower].

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

(Contd. on page 2)

Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)

P310 Immediately call a POISON CENTER/doctor.
P321 Specific treatment (see on this label).
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

## 3 Composition/information on ingredients

- · Mixtures
- · **Description:** Mixture of substances listed below with nonhazardous additions.
- · Dangerous components: Void
- · Additional information: For the wording of the listed hazard phrases refer to section 16.

## 4 First aid measures

- · Description of first aid measures
- · After inhalation:

Move to fresh air. If breathing is difficult, give oxygen. Perform artificial respiration if breathing has stopped. Get medical attention immediately. Keep patient warm

- · After skin contact: Immediately flush skin with plenty of cool soapy water. Get medical attention
- · After eye contact:

Immediately flush eyes with plenty of water for several minutes, while holding eyelids apart. Get medical attention immediately.

- · After swallowing: If ingestion occurs, flush mouth with water, and seek medical attention immediately.
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed Seek medical advice after large-scale exposure.

#### 5 Firefighting measures

- · Extinguishing media Water spray, foam, carbon dioxide, or dry chemical.
- · Suitable extinguishing agents: None.
- · Special hazards arising from the substance or mixture

Fumes of carbon dioxide, nitrogen oxides, and carbon monoxide may be evolved in a fire.

· Advice for firefighters: Usual. Avoid inhalation of smoke; use a respirator and protective suit.

#### 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions: None.
- · Methods and material for containment and cleaning up:

Wipe up spills with rags to prevent material from being dispersed, and discard in trash.

Dispose of container and unused contencts in accordance with federal, state, and local requirements.

· Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)

(Contd. of page 2)

#### 7 Handling and storage

- · Precautions for safe handling:
- · Handling:

Gel is slightly acidic. Wear suitable gloves and protective glasses during use. Keep out of eyes or contact with mucosal tissues. After incidental contact, wash and rinse with water. Endo-Eze MTAFlow Gel must be kept well-sealed, and free of contamination by pyrogens. Only trained personnel may use this product exclusively for dental purposes.

· Conditions for safe storage, including any incompatibilities

Store at room temperature (25C/77F). Keep bottle tightly closed. Inappropriate storage will shorten the shelf life.

- · Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

#### 8 Exposure controls/personal protection

- · Control parameters
- · Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Exposure controls
- · Individual protection measures, such as personal protective equipment
- · General protective and hygienic measures:

Use in a well-ventilated area.

Wear gloves as per universal precautions in dentistry. Wash your hands after use. No eating, drinking, or smoking while handling this product. Wash clothing with spills of material on clothing.

· Respiratory protection:

Avoid inhalation

Dust mask

Respirator

- · Hand protection Wear latex or nitrile gloves
- · Eye/face protection Wear safety glasses or face shield.

## 9 Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information

Physical state
Colour:
Pale Yellow
Odour:
Odourless
Odour threshold:
Not determined.

• Melting point/freezing point: 0 °C • Boiling point or initial boiling point and boiling range 100 °C

· Flammability Not applicable.

· Lower and upper explosion limit

· **Lower:** Not determined.

(Contd. on page 4)

Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)

	(Contd. of pag
Upper:	Not determined.
Flash point:	Not applicable.
Decomposition temperature:	Not determined.
pH at 20 °C	3-4
Viscosity:	
Kinematic viscosity	Not determined.
Dynamic:	Not determined.
Solubility	
water:	Fully miscible.
Partition coefficient n-octanol/water (log value)	Not determined.
Vapour pressure:	Not determined.
Density and/or relative density	
Density at 20 °C:	$1.1 \text{ g/cm}^3$
Relative density	Not determined.
Vapour density	Not determined.
Other information	
Appearance:	
Form:	Gel
Important information on protection of health a environment, and on safety. Auto-ignition temperature: Explosive properties:	Product is not selfigniting. Product does not present an explosion hazard.
Change in condition	Trouder does not present an expression nazara.
Evaporation rate	Not determined.
Information with regard to physical hazard classes	
Explosives	Void
Flammable gases	Void
Aerosols	Void
Oxidising gases	Void
Gases under pressure	Void
Flammable liquids	Void
Flammable solids	Void
Self-reactive substances and mixtures	Void
Pyrophoric liquids	Void
Pyrophoric solids	Void
Self-heating substances and mixtures	Void
Substances and mixtures, which emit flammable gas	
in contact with water	Void
Oxidising liquids	Void
Oxidising solids	Void
Organic peroxides	Void
Corrosive to metals	Void
Desensitised explosives	Void

## 10 Stability and reactivity

- · Reactivity No further relevant information available.
- · Chemical stability Stable under recommended conditions.
- · Possibility of hazardous reactions: No dangerous reactions known.
- · Conditions to avoid: Heat
- · Incompatible materials: Do not mix with material other than MTA Flow Powders.

(Contd. on page 5)

Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)

(Contd. of page 4)

· Hazardous decomposition products:

Decomposition products of the minor components of gel include carbon dioxide, nitrogen oxiders, and carbon monoxide.

## 11 Toxicological information

- · Information on hazard classes as defined in Regulation (EC) No 1272/2008
- · Acute toxicity Based on available data, the classification criteria are not met.
- · Skin corrosion/irritation None.
- · Serious eye damage/irritation None.
- Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- · Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- · Reproductive toxicity Based on available data, the classification criteria are not met.
- · STOT-single exposure Based on available data, the classification criteria are not met.
- · STOT-repeated exposure Based on available data, the classification criteria are not met.
- · Aspiration hazard Based on available data, the classification criteria are not met.
- · Additional toxicological information: None.
- Information on other hazards
- · Endocrine disrupting properties

None of the ingredients is listed.

## 12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.
- Endocrine disrupting properties For information on endocrine disrupting properties see section 11.
- · Other adverse effects

Toxicity, Environmental Fate, Physical/Chemical Data, or Other Data Supporting Environmental Hazard Statements; Discard in Trash

## 13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Dispose of contents/container in accordance with international, federal, state, and local regulations.

- · Uncleaned packaging:
- · Recommendation: None.

GB

Printing date 23.02.2024 Version number 1 Revision: 23.02.2024

Trade name: Endo-Eze<sup>TM</sup> MTAFlow<sup>TM</sup> Gel (Gray and White)

(Contd. of page 5)

Transport information	
UN number or ID number ADR, IMDG, IATA	not regulated
UN proper shipping name ADR, ADN, IMDG, IATA	not regulated
Transport hazard class(es)	
ADR, ADN, IMDG, IATA Class	not regulated
Packing group ADR, IMDG, IATA	not regulated
Environmental hazards:	Not applicable.
Special precautions for user	Not Applicable
Maritime transport in bulk according instruments	<b>g to IMO</b> Not applicable.
UN "Model Regulation":	not regulated

## 15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- · Chemical safety assessment: A chemical safety assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing SDS: Environmental, Health, and Safety
- · Contact: Customer Service
- · Abbreviations and acronyms:

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Skin Corr. 1A: Skin corrosion/irritation – Category 1A

\* \* Data compared to the previous version altered.