



## Safety data sheet according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

### 1 Identification of the substance/mixture and of the company/undertaking

- **Product identifier**
- **Trade name:** Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good
- **Article number:**  
SDS 505-1000.03R01, 1011957, 1012434, 5822-CA, 5822-US, 6006-CA, 6006-US, 6006-U, 1014114, 1014116, 15441, 5822A-US, 6004A-US, 6006A-US, 6008A-US, S6004A-US
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** At-Home Professional-Grade Teeth Whitening Gel
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ultradent Products Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com  
(800) 552-5512
- **EC Responsible Person**  
Ultradent Products GmbH  
Am Westhover Berg 30  
51149 Cologne Germany  
Email: infoDE@ultradent.com  
Office Phone: +49(0)2203-35-92-0
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**  
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

### 2 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**  
The product is not classified, according to the GB CLP regulation.
- **Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void

### 3 Composition/information on ingredients

- **Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

(Contd. on page 2)

# Safety data sheet

## according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 1)

· <b>Dangerous components:</b>		
CAS: 56-81-5 EINECS: 200-289-5	Glycerin substance with a Community workplace exposure limit	1-10%
CAS: 7722-84-1 EINECS: 231-765-0	Hydrogen Peroxide ⚠ Ox. Liq. 1, H271; ⚠ Skin Corr. 1A, H314; ⚠ Acute Tox. 4, H302; Acute Tox. 4, H332 Specific concentration limits: Ox. Liq. 1; H271: $C \geq 70\%$ Ox. Liq. 2; H272: $50\% \leq C < 70\%$ Skin Corr. 1A; H314: $C \geq 70\%$ Skin Corr. 1B; H314: $50\% \leq C < 70\%$ Skin Irrit. 2; H315: $35\% \leq C < 50\%$ Eye Dam. 1; H318: $C \geq 8\%$ Eye Irrit. 2; H319: $5\% \leq C < 8\%$ STOT SE 3; H335: $C \geq 35\%$	$\geq 1$ -<5%
CAS: 8006-90-4 EINECS: 282-015-4	Oils, Peppermint ⚠ Asp. Tox. 1, H304; ⚠ Aquatic Chronic 2, H411; ⚠ Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317	$\geq 0.1$ -<0.25%

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

### 4 First aid measures

- **Description of first aid measures**
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters:**
- **Protective equipment:** No special measures required.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

GB

(Contd. on page 3)

# Safety data sheet

## according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 2)

### 7 Handling and storage

- **Precautions for safe handling:** See product labeling.
- **Information about fire - and explosion protection:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Control parameters**

- **Ingredients with limit values that require monitoring at the workplace:**

#### 56-81-5 Glycerin

WEL Long-term value: 10 mg/m<sup>3</sup>

#### 7722-84-1 Hydrogen Peroxide

WEL Short-term value: 2.8 mg/m<sup>3</sup>, 2 ppm

Long-term value: 1.4 mg/m<sup>3</sup>, 1 ppm

- **Additional information:** The lists valid during the making were used as basis.
- **Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:** Wash hands before breaks and at the end of work.
- **Respiratory protection:** Not required
- **Hand protection**  
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**  
The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**  
The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye/face protection** Not required
- **Body protection:** Protective work clothing

### 9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Physical state** Solid
- **Colour:** White Opaque
- **Odour:** Mint

(Contd. on page 4)

## Safety data sheet according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 3)

· <b>Odour threshold:</b>	<i>Not determined</i>
· <b>Melting point/freezing point:</b>	<i>Undetermined</i>
· <b>Boiling point or initial boiling point and boiling range</b>	<i>Undetermined</i>
· <b>Flammability</b>	<i>Not determined</i>
· <b>Lower and upper explosion limit</b>	
· <b>Lower:</b>	<i>Not determined</i>
· <b>Upper:</b>	<i>Not determined</i>
· <b>Flash point:</b>	<i>Not applicable</i>
· <b>Decomposition temperature:</b>	<i>Not determined</i>
· <b>pH at 20 °C</b>	<i>5-7</i>
· <b>Viscosity:</b>	
· <b>Kinematic viscosity</b>	<i>Not applicable</i>
· <b>Dynamic:</b>	<i>Not applicable</i>
· <b>Solubility</b>	
· <b>water:</b>	<i>Insoluble.</i>
· <b>Partition coefficient n-octanol/water (log value)</b>	<i>Not determined</i>
· <b>Vapour pressure:</b>	<i>Not applicable</i>
· <b>Density and/or relative density</b>	
· <b>Density at 20 °C:</b>	<i>1.3 g/cm<sup>3</sup></i>
· <b>Relative density</b>	<i>Not determined</i>
· <b>Vapour density</b>	<i>Not applicable</i>

· <b>Other information</b>	
· <b>Appearance:</b>	
· <b>Form:</b>	<i>Gel</i>
· <b>Important information on protection of health and environment, and on safety.</b>	
· <b>Ignition temperature:</b>	<i>Product is not selfigniting.</i>
· <b>Explosive properties:</b>	<i>Product does not present an explosion hazard.</i>
· <b>Change in condition</b>	
· <b>Evaporation rate</b>	<i>Not applicable.</i>

· <b>Information with regard to physical hazard classes</b>	
· <b>Explosives</b>	<i>Void</i>
· <b>Flammable gases</b>	<i>Void</i>
· <b>Aerosols</b>	<i>Void</i>
· <b>Oxidising gases</b>	<i>Void</i>
· <b>Gases under pressure</b>	<i>Void</i>
· <b>Flammable liquids</b>	<i>Void</i>
· <b>Flammable solids</b>	<i>Void</i>
· <b>Self-reactive substances and mixtures</b>	<i>Void</i>
· <b>Pyrophoric liquids</b>	<i>Void</i>
· <b>Pyrophoric solids</b>	<i>Void</i>
· <b>Self-heating substances and mixtures</b>	<i>Void</i>
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	<i>Void</i>
· <b>Oxidising liquids</b>	<i>Void</i>
· <b>Oxidising solids</b>	<i>Void</i>
· <b>Organic peroxides</b>	<i>Void</i>
· <b>Corrosive to metals</b>	<i>Void</i>
· <b>Desensitised explosives</b>	<i>Void</i>

GB

(Contd. on page 5)

## Safety data sheet according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 4)

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

### 11 Toxicological information

- **Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

#### ATE (Acute Toxicity Estimates)

Oral	LD50	19,335 mg/kg
Inhalative	LC50/4 h	>382 mg/l

#### 56-81-5 Glycerin

Oral	LD50	7,750 mg/kg (guinea pig)
		4,100 mg/kg (mouse)
		5,570 mg/kg (rat)
		27,000 mg/kg (rabbit)
Dermal	LC50 Fish	>5,000 mg/l (Fish)
	LD50	>21,900 mg/kg (rat)
Inhalative	LC50/4 h	10,000 mg/kg (rabbit)
		>0.1425 mg/l (rat)

#### 7722-84-1 Hydrogen Peroxide

Oral	LC50 Fish	16.4 mg/l (Fish)
------	-----------	------------------

#### 8006-90-4 Oils, Peppermint

Oral	LD50	2,490 mg/kg (mouse)
		2,426 mg/kg (rat)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

(Contd. on page 6)

# Safety data sheet

## according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 5)

· **Information on other hazards**

· **Endocrine disrupting properties**

7681-49-4	Sodium Fluoride	List II
-----------	-----------------	---------

### 12 Ecological information

· **Toxicity**

· **Aquatic toxicity:**

**56-81-5 Glycerin**

EC50	>10,000 mg/kg (Bacteria)
------	--------------------------

**7722-84-1 Hydrogen Peroxide**

EC50	1.38 mg/l (Algae)
------	-------------------

	2.4 mg/l (daphnia)
--	--------------------

· **Persistence and degradability** No further relevant information available.

· **Bioaccumulative potential** No further relevant information available.

· **Mobility in soil** No further relevant information available.

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable

· **vPvB:** Not applicable

· **Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.

· **Other adverse effects**

· **Additional ecological information:**

· **General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

### 13 Disposal considerations

· **Waste treatment methods**

· **Recommendation**

Dispose of contents/container in accordance with international, federal, state, and local regulations.

· **Uncleaned packaging:**

· **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

· **UN number or ID number**

· **ADR, IMDG, IATA** not regulated

· **UN proper shipping name**

· **ADR, ADN, IMDG, IATA** not regulated

· **Transport hazard class(es)**

· **ADR, ADN, IMDG, IATA**

· **Class** not regulated

· **Packing group**

· **ADR, IMDG, IATA** not regulated

(Contd. on page 7)

## Safety data sheet according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 6)

- |  |                |
|--|----------------|
| · <b>Environmental hazards:</b>                                  | Not applicable |
| · <b>Special precautions for user</b>                            | Not Applicable |
| · <b>Maritime transport in bulk according to IMO instruments</b> | Not applicable |
| · <b>UN "Model Regulation":</b>                                  | not regulated  |

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

- **Poisons Act**

- **Regulated explosives precursors**

7722-84-1	Hydrogen Peroxide	12%
-----------	-------------------	-----

- **Regulated poisons**

None of the ingredients is listed.

- **Reportable explosives precursors**

7757-79-1	Potassium Nitrate	Listed
-----------	-------------------	--------

- **Reportable poisons**

1310-58-3	Potassium Hydroxide	17% of total caustic alkalinity
-----------	---------------------	---------------------------------

7681-49-4	Sodium Fluoride	Listed
-----------	-----------------	--------

- **Directive 2012/18/EU**

- **Named dangerous substances - ANNEX I** None of the ingredients is listed.

- **Chemical safety assessment:** A chemical safety assessment has not been carried out.

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases from Section 3**

H271 May cause fire or explosion; strong oxidiser.

H272 May intensify fire; oxidiser.

H302 Harmful if swallowed.

H304 May be fatal if swallowed and enters airways.

H314 Causes severe skin burns and eye damage.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H335 May cause respiratory irritation.

H411 Toxic to aquatic life with long lasting effects.

- **Department issuing SDS:** Regulatory Compliance

- **Contact:** Customer Service

(Contd. on page 8)

## Safety data sheet according to UK REACH

Printing date 28.05.2026

Version number 1

Revision: 28.05.2026

**Trade name: Opalescence™ Go™ Advanced Tooth Whitening 15% - Finished Good**

(Contd. of page 7)

· **Abbreviations and acronyms:**

*ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*IATA: International Air Transport Association*

*GHS: Globally Harmonised System of Classification and Labelling of Chemicals*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*LC50: Lethal concentration, 50 percent*

*LD50: Lethal dose, 50 percent*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*NIOSH: National Institute for Occupational Safety*

*ATE: Acute toxicity estimate values*

*Ox. Liq. 1: Oxidizing liquids – Category 1*

*Acute Tox. 4: Acute toxicity – Category 4*

*Skin Corr. 1A: Skin corrosion/irritation – Category 1A*

*Skin Irrit. 2: Skin corrosion/irritation – Category 2*

*Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

*Skin Sens. 1: Skin sensitisation – Category 1*

*Asp. Tox. 1: Aspiration hazard – Category 1*

*Aquatic Chronic 2: Hazardous to the aquatic environment - long-term aquatic hazard – Category 2*

· **\* Data compared to the previous version altered.**

GB